



iMedColorRx

Color pharmacy labels printed on demand

The iMedColorRx Advantage



Alerts & Special Handling

Print color alerts and special handling indicators directly on each label, preventing labeling errors.



Streamlined Pharmacy Operations

Reduce manual intervention and cross-check procedures for the manual application of alert labels.



Pediatric/Adult Color-Coding

Reduce medication administration errors by color-coding pediatric and adult doses.



Tall Man Lettering

Reduce medication errors for look-alike/sound-alike drugs through user-defined Tall Man lettering.



Dynamic Printing

Reduce costs by eliminating separate alert labels and associated inventory management.



Output Management

Leverage cut-to-size and “roll-up” logic to eliminate unneeded white space and label waste.

TAYLOR
HEALTHCARE

7SXX 6702

STAT **PROTECT FROM LIGHT** **DO NOT REFRIGERATE**

Acct: 99999999 DOB: 99/99/9999 5Y F

PatientName, First

ORD#: 9999 BAG: 999 DUE: XX XX IV

Norepinephrine 4mg = 4mL - DEXTROSE 5% - Water 250 ML

Rate: CONTINUOUS Freq: CONTINUOUS

Initiate at .1 MCG/MIN. Adjust bt .5 MCG/MIN Q3 MIN to Maintain SBP greater than 90. Usual range .2-7.1 MCG/MIN. Do NOT begin after 23:59 on 6/11/13

Generic Name Addt 01	Mod Dose Addt 01
Generic Name Addt 02	Mod Dose Addt 02
Generic Name Addt 03	Mod Dose Addt 03
Generic Name Addt 04	Mod Dose Addt 04
Generic Name Addt 05	Mod Dose Addt 05

Easily differentiate between pediatric and adult patients

Eliminate manual auxiliary label processes for special handling conditions

Enhance communications with colorfast, smudge-proof labels

Reinforce alerts and warnings (e.g., STAT, high risk, chemo)

User-defined Tall Man table distinguishes between look-alike/sound-alike drugs

Highlight additional warnings

Highlight critical information

TAYLOR
HEALTHCARE

7SXX 6702

STAT **CHMO** **HAZARDOUS** **HIGH ALERT**

Acct: 99999999 DOB: 99/99/9999 27Y F

PatientName, First

ORD#: 9999 SYR: 999 DUE: XX XX *START* SYRINGE

vinBLASTine 15.3 MG = 15.3 ML TOTAL VOLUME

Route: IT Freq: ON CALL

Do NOT begin after 23:59 on 99/99/99

WEAR GLOVES
FATAL IF GIVEN INTRATHECALLY VESSICANT

Generic Name Addt 01	Mod Dose Addt 01
Generic Name Addt 02	Mod Dose Addt 02
Generic Name Addt 03	Mod Dose Addt 03
Generic Name Addt 04	Mod Dose Addt 04

Barcode: [Barcode] Printed: 99/99/9999 99.99 RPH: Tech: _____

How It Works



Pharmacy HIS

Pharmacy HIS sends data feed or HL7 RDS messages to Taylor Healthcare's EPS server for processing pharmacy labels.

→
Data/
HL7



EPS Server

EPS server software transforms the data, processes color and routing logic, renders output, and securely sends to color inkjet label printers.

→
TCP/
IP



Color Printer

Customer prints crisp, vibrant medication, IV and piggyback labels that won't run or smudge. Internal cutter allows optional cut-to-size functionality.